

Pressureplate not working on EVO (LED states)

Please see below in green the standard behaviour of the LEDs in the white frame on top of the bottom plate.

LED 1 - Green/Red (not relevant for EVO)		LED 2 - Blue (Data)		LED 3 - Yellow (Power)	
On	Old firmware version (first 10 EVOs)	On	USB connected and Windows driver detects the plate	5 sec active on startup	Plate has power and starts.
		Pulse	When software is in measure mode and sensor detects weight	Pulse	Power on.

Troubleshooting

LED	is	but should be	Solution
Yellow (LED 3)	Off	On	<ul style="list-style-type: none">• Check Power Button on plate. The power button is inside the housing and has a quick access with just a few screws on top.• Check USB cabling between computer and pressure plate
Blue (LED 2)	Off	On or Pulse	
RED (LED 1)	Blinking		Firmware before April 2020 If the pressure plate is not working, exchange it with a new firmware one. (If the pressure plate works, don't exchange anything)